

Project	
Context Sensitive Help	F1
Find in Project	F3
New	Ctrl+N
Open	Ctrl+O
Save Project	Ctrl+S
Print	Ctrl+P
Switch Windows in Current Folder	Ctrl+Tab
Return to Last Window	Ctrl+Alt+Z
Toggle Project Explorer View	Alt+1
Toggle Command Window View	Alt+3
Toggle Script Window View	Alt+Shift+3
Activate Code Builder	Alt+4
Toggle Messages Log View	Alt+6
Toggle Smart Hint Log	Alt+7
Toggle Object Manager	Alt+8
Toggle Apps Gallery	Alt+9

Scale and Pan	
Rescale to Show All Data	Ctrl+R
Scale (Up=In, Down=Out) the X axis	Z+Wheel
Scale (Up=In, Down=Out) the Y axis	Shift+Z+Wheel
Panning	Z+Mouse
Panning along the X axis	X+Mouse
Panning along the Y axis	Shift+X+Mouse

Tools	
Theme Organizer	F7
Fitting Function Builder	F8
X-Function Builder	F10
Origin Central	F11
Import Wizard	Ctrl+3
Customise Toolbars	Ctrl+T
Options	Ctrl+U
Open NLFit	Ctrl+Y
Open XF Dialogs with Theme Loaded	Shift+menu with popup <theme>

Column	
Set Column Values	Ctrl+Q
Set Multiple Column Values	Ctrl+Shift+Q
Copy Column Formula	Ctrl+J
Copy All (Formula, Label Rows, Values)	Ctrl+Shift+A

Editing	
Toggle In-Place Edit Mode	F2
Clear	Delete
Select All	Ctrl+A
Copy	Ctrl+C
Copy Full Precision	Ctrl+Alt+C
Copy then Delete	Ctrl+X
Paste	Ctrl+V
Paste Link	Ctrl+Alt+V
Undo	Ctrl+Z
Find	Ctrl+F
Replace	Ctrl+H
Toggle Button Edit Mode	Ctrl+Alt+B

Worksheet	
Worksheet Properties	F4
Copy Including Labels	Ctrl+Shift+C
Add New Columns	Ctrl+D
Go To	Ctrl+G
Go to First/Last Cell in Worksheet	Ctrl+ Home/End
Switch Worksheets in Workbook	Ctrl+PageUp/ PageDown

Matrix	
Set Matrix Value	Ctrl+Q
Row/Column Index View	Ctrl+Shift+C
Data View	Ctrl+Shift+D
Image View	Ctrl+Shift+I
Row/Column Value View	Ctrl+Shift+X
Toggle Thumbnail View	Ctrl+Shift+N
Toggle ROI Tools	Ctrl+Shift+T

Import and Export	
Import ASCII	Ctrl+K
Import Wizard	Ctrl+3
Reimport Directly	Ctrl+4
Export Graph	Ctrl+G

Zoom	
Zoom (Up =In, Down = Out)	Ctrl+Wheel
100% View	Ctrl+W
Step Zoom In Tool	Ctrl+I
Step Zoom Out	Ctrl+M



**ORIGIN® 2017**  
Graphing & Analysis

Graph	
Open Plot Details	F2
Single Point customization	Ctrl+ Double-click
Select Layer	Alt+Click
Data Reader	3
Annotation	***
In Annotation Mode: Change Annotation Syntax	Tab
Data Selector	Ctrl+4
In Data Selection Mode: Toggle Selection on Active Plot/All Plots	Tab
In Data Selection Mode: Change Select Mode	Space
Draw Data	Ctrl+5
Text Tool	Ctrl+6
Insert Variable (Label)	Ctrl+H
Copy Page	Ctrl+J
Update Layer Legend	Ctrl+L
Symbol Map (Label)	Ctrl+M
Rescale Graph Layer	Ctrl+R
Full Screen View	Ctrl+Shift+J

**OriginLab®**

[www.originlab.com](http://www.originlab.com)

*25+ years serving the scientific and engineering community.*

# Origin 2017 Toolbar Quick-Reference

**Standard Tools**

- Create function plots
- Zoom in and out with axis rescaling
- Read, annotate, and measure data points
- Select a region of data on one, or all plots for analysis
- Mask and unmask data points on one, or all plots
- Add annotation
- Add text, arrows, lines, and shapes
- Insert an equation, MS Word, MS Excel, Floating Graph or Worksheet
- Rotate 3D graphs

**Graph**

- Save a template
- Import wizard
- Recalculate analysis results
- Duplicate graph with new worksheets
- Perform batch processing
- Send graphs to PowerPoint
- Create video from graphs
- Rescale graph to show all data
- Rescale the X Axis
- Rescale the Y Axis
- Merge cells
- Duplicate graph with new sheet/book

**Style**

- Color and line style options

**Worksheet Data**

- Data Filter

**Format**

- Font and text formatting options

**2D Graph Types**

- Add color scale from color map
- Create a bubble scale
- Create or update legend
- Multi-axis, multi-panel, waterfall
- Radar, polar, vector, ternary, windrose

**3D Graph Types**

- 3D Bars
- 3D Surface
- 3D Scatters
- Contours

**OriginLab**  
www.originlab.com

25+ years serving the scientific and engineering community.